





PharmaMV PAT, Advanced Control and Optimization Platform



Andy Mitchell
Product Development Director
APMF 2015



Advanced Process Modelling Forum

Agenda

- **Introduction to Perceptive Engineering**
- **PharmaMV Platform**
- **Perceptive Capabilities**
- **Look ahead to use of PharmaMV & gCRYSTAL models**

- **Spin-out 2003 Manchester University**
 - Department of Advanced Control
- **Engineers practical APC experience**
 - Average 8 years practical experience
 - Principal Team Leaders >15 years experience
- **Offices in UK and Singapore**
- **Academic alliances**
 - Universities of Cambridge, Manchester, Newcastle, TU Delft, UCL
 - CMAC & ICES – Alliance of several academic and industrial partners in continuous crystallization
 - CPI – Centre for Process Innovation



- **What We Do:**

- Perceptive Engineering develop and use a software toolset combined with engineering consultancy to create systems that optimise industrial manufacturing processes.
- Introduction of control and monitoring technologies and techniques into a commercial offering for most non oil and gas manufacturing sectors.

- **Our Approach to R&D**

- Translation of academic research into market driven and robust industrial solutions.
- Collaborative R&D leverages resources and skillsets from multiple parties, each deriving benefits, to meet their own distinctive needs.



Nestlé



Wyeth



Lilly



Perceptive Engineering Software Family



Off-line development, analysis, visualisation and assessment:
find out how capable your plant can be



On-line process improvement:
get the most from your process, all the time




Industry-specific enhanced toolsets:
designed to meet unique challenges

The PharmaMV Platform Philosophy

- **A software platform for Advanced Process Control in the Pharmaceutical Industry:**
- **Some of the challenges....**
 - Management of “Instrument and Data” Integrity
 - Capability to deal with Batch and Continuous Processes
 - Interconnectivity – Use of existing industrial standards.
 - Delivering Fault tolerant Monitoring, Control, Optimization.
 - Traceable User Actions and 21CFR Part 11 Compliance.
- **Blending the Pharmaceutical Scientist and Automation Engineering disciplines**

Manufacturing Intelligence Software for Real Time Analysis and Control



Data Import/Export



Data Quality Monitoring



Data Analysis



SPC Monitoring



Multivariate Modelling



Multivariate Process Monitoring




Multivariate MPC




**Continuous / Batch
Optimisation**

Before we go any further...

Setting the scene on the word "Model"




- **Stochastic Models for**
 - Dynamic Control
 - Predictions
 - Calibration etc.
- **Created from**
 - Designed Plant Tests
 - Historical Process Data




PharmaMV Platform Architecture

Real-Time Connectivity




PharmaMV Platform

Real-Time Manager and Spectral Data




Univariate Process Monitoring Continuous Process Verification

- **Real-time abnormal operation detection**
 - Statistical rules
 - Control Limits update based on Product / Lot / Operation of process




Multivariate Process Monitoring Continuous Process Verification



Diagnose Process Fault

Process Abnormal Operation



Guided Analysis for both
Operator and Engineer.

Automated Plant Testing

Design of Experiments / Process Response Testing

Modelling Facilities

Specification | Coefficients | Evaluation

Details | Batch Details | Models | Record

Predictor (X)... Response (Y)... Options... Save... Save Record Set... Load... New

Block 1000

Design of Experiments

Response (Y) Signal	Minimum	Maximum
1.ME Output 1	50	60
2.ME Output 2	40	70
3.ME Output 3	50	75
4.ME Output 4	0	10

Predictor (X) Signal

Predictor (X) Signal	Minimum
1.AC Input 1	50
2.AC Input 2	50
3.AC Input 3	0
4.AC Input 4	0

Configuration for 1000

Design of Experiments

Number of Points: 32

Time to Steady State: 10.0s

Steady State Multiples of Point: 2

Step Steady State Fractions: 0.5

Step in DoE Direction Only: No

DoE Point	Input 1	Input 2	Input 3	Input 4
1	50	40	50	2.5
2	50	40	50	5
3	50	40	60	2.5
4	50	40	60	5
5	50	50	50	2.5
6	50	50	50	5
7	50	50	60	2.5
8	50	50	60	5
9	55	50	60	5
10	55	50	60	5
11	55	50	75	5
12	55	50	75	5
13	55	70	50	5
14	55	70	50	5

Expected Relationships



	Input 1	Input 2
Output 1	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Output 2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Output 3	<input type="checkbox"/>	<input type="checkbox"/>
Output 4	<input type="checkbox"/>	<input type="checkbox"/>

Define the boundaries


Setup Test Conditions

Perform and Monitor Tests

Process Analytical Technology (PAT) Data Processing / Analysis and Calibration Models



- **Using Process Analytical Technology (PAT) in control offers great possibilities**
 - Using modern PAT devices capable of 1-2 second or sub-second measurement rates, real-time control based on PAT measurement is a reality
 - Sensor calibration models can give real time inference of product property



PharmaMV Platform


Pharma and CPG Applications

- **APC for API Manufacture**
 - Crystallisation
 - Continuous and Batch
 - Solvent Extraction
- **APC for Solid Dose Form**
 - Blending, Granulation, Drying, Coating and Tableting
- **APC for Fermentation**
 - Run to Run Optimisation
- **APC for Drying**
 - Nutritional Powders



PharmaMV Platform


Batch Crystallisation Concentration Prediction



PharmaMV Platform


Batch Crystallisation APC System

Batch Crystallization Real-Time APC System



An integrated platform for improved cool-down crystallization

- Automated Meta-stable zone identification for product development and scale-up.
- Validated Real-time concentration prediction.
- Robust super-saturation control.



PharmaMV Platform


High Shear Wet Granulation (HSWG)

- **Monitoring and/or control of High-Shear Wet Granulation**
 - NIR PAT modelled to generate latent variable score trajectories in real-time
 - Batch unfolded techniques to control and/or monitor the process based on these trajectories
- **Objectives:**
 - Control the process to a “golden batch” trajectory
 - Respond to variations in raw material and other factors

PharmaMV Platform

High Shear Wet Granulation (HSWG)

HSWG Real-Time APC System




An integrated platform for HSWG control and monitoring



- The latent variable scores describe the key features of the NIR spectrum.
- A batch-end point controller is used to adjust the batch trajectory at key decision points. The controller automatically switches models based on the batch phase.

PharmaMV Platform

High Shear Wet Granulation (HSWG) for Operators




Continuous Crystallisation Continuous Oscillatory Baffle Reactor



• Main features:

- Data acquisition, fail-safe automation interface and operator advisory display for
 - Thermocouples
 - Oscillator
 - Pumps
- Data acquisition and interface to PAT devices including
 - FBRM
 - FTIR
 - Raman
 - React IR
- Vessel Temperature Control display




PharmaMV Platform

Continuous Crystallisation (CSTR/MSMPR/Anti Solvent)

- PharmaMV Platform

- Flexible design to allow unit to be operated in different modes dependant on requirements



Nutritional Powders Solution APC & Optimisation

**Data Quality
Monitor**

KPIs


**Thermodynamic
Calculations**

Soft sensors

**Evaporator/ Dryer
Start-up**

**APC &
Optimisation**

REAL TIME



perceptiveAPC

Controlled Variables	Target	SP	Process Value	Set Point	Control Value	High Limit	Low Limit
1. Inlet Air Humidity	45.00	45.00	45.00	45.00	45.00	45.00	45.00
2. Inlet Air Temperature	45.00	45.00	45.00	45.00	45.00	45.00	45.00
3. Inlet Air Flow	45.00	45.00	45.00	45.00	45.00	45.00	45.00
4. Inlet Air Pressure	45.00	45.00	45.00	45.00	45.00	45.00	45.00
5. Inlet Air Humidity	45.00	45.00	45.00	45.00	45.00	45.00	45.00
6. Inlet Air Temperature	45.00	45.00	45.00	45.00	45.00	45.00	45.00
7. Inlet Air Flow	45.00	45.00	45.00	45.00	45.00	45.00	45.00
8. Inlet Air Pressure	45.00	45.00	45.00	45.00	45.00	45.00	45.00
9. Inlet Air Humidity	45.00	45.00	45.00	45.00	45.00	45.00	45.00
10. Inlet Air Temperature	45.00	45.00	45.00	45.00	45.00	45.00	45.00
11. Inlet Air Flow	45.00	45.00	45.00	45.00	45.00	45.00	45.00
12. Inlet Air Pressure	45.00	45.00	45.00	45.00	45.00	45.00	45.00
13. Inlet Air Humidity	45.00	45.00	45.00	45.00	45.00	45.00	45.00
14. Inlet Air Temperature	45.00	45.00	45.00	45.00	45.00	45.00	45.00
15. Inlet Air Flow	45.00	45.00	45.00	45.00	45.00	45.00	45.00
16. Inlet Air Pressure	45.00	45.00	45.00	45.00	45.00	45.00	45.00
17. Inlet Air Humidity	45.00	45.00	45.00	45.00	45.00	45.00	45.00
18. Inlet Air Temperature	45.00	45.00	45.00	45.00	45.00	45.00	45.00
19. Inlet Air Flow	45.00	45.00	45.00	45.00	45.00	45.00	45.00
20. Inlet Air Pressure	45.00	45.00	45.00	45.00	45.00	45.00	45.00




• BENEFITS

- Stabilise process & reduce variability
- Improve quality compliance
- Increase yield and capacity
- Reduce specific energy consumption




Nutritional Powders Solution APC & Optimisation

- The Foundation of our Solution




Nutritional Powders Solution Evaporator and Dryer Start-up

Water → Product start-up with APC



Multiples of 16s (1 hour span)



Multiples of 32s (1 hour span)



Process Monitoring, Control and Optimisation

Look ahead to new approaches

- **Process Monitoring / Control / Simulation and Optimisation**
 - Current Models are data driven
 - Involves plant testing
 - Cost of raw materials
 - Availability of plant for testing
 - Existing historical data could be in a narrow band of operation
 - Could
 - Plant testing be performed on a mechanistic model?
 - A process be monitored / controlled directly via a mechanistic model?
 - Increase the use of multivariable control / monitoring at design time?
 - Process capability increase at design time?

PharmaMV with gCRYSTAL Process Units

- **PharmaMV using gCRYSTAL as a process unit**
 - Process testing and advanced control design




- **Potential Application & Benefits**
 - Selection of control parameters, initial process testing and model design on a model
 - Reduction on material required to build models
 - Reduction on time required on site, process availability
 - Reduced control/monitor commissioning time.

PharmaMV with gCRYSTAL

Mixed / Hybrid approach

- **Mixed / Hybrid approach to control and monitoring.**



- **Potential Application & Benefits**
 - Reduced or minimal process testing required
 - Models available for variables that may not have been possible before

Any Questions... Just Ask!



Andy Mitchell
Product Development Director
amitchell@perceptiveapc.com