



Improving Process Understanding through Crystallization Modeling

Shane T. Grosser¹

Eric Sirota¹, Jochen Schoell² and Lorenzo Codan²

2013 Advanced Process Modeling Forum

New York, NY


June 5th, 2013

¹Chemical Process Development and Commercialization, Merck, Rahway, NJ,

²Chemical Process Development and Commercialization, Merck, Wertheim, Switzerland

- ❑ Crystallization is frequently employed to produce solid material fit for formulation from chemical synthesis steps
- ❑ Oral absorption and formulation performance can be closely tied to particle physical attributes (i.e. particle size and shape)
 - Historically dry milling has been employed to reduce particle size to improve adsorption performance. Significant cost and IH concerns driving a desire for direct crystallization with minimal dry milling required.
- ❑ Overall process time cycle and yearly capacity can also be negatively impacted by slowly filtering compounds
 - Filtration performance is frequently correlated to particle size and shape and frequently benefits from larger, more growth dominated crystallization processes.
- ❑ Crystallization design plays a key role in balancing many of these particle attribute considerations

- ❑ Merck maintains a dedicated crystallization laboratory to support small molecule crystallization design and development
- ❑ Staffed by 9 engineers, the X-lab focuses on support for:
 - Crystallization development support of process design
 - Drug substance processability assessment
 - ❑ Solubility screening
 - ❑ Wet-milling performance evaluation
 - ❑ Evaluation of filtration characteristics
 - ❑ Sensitivity to particle breakage
 - Intermediate and drug substance polymorph screening
 - Supply stabilization support and process improvement
 - Crystallization technology development
- ❑ X-lab operates as a center of excellence with significant subject matter expertise and advanced development and analytical tools.



☐ Particle Characterization

- SEM
- Optical Microscope
- Sympatec
- Microtrac
- Surface Area
- Powder density

☐ Wet Milling

- Rotor stator milling
- Media milling

☐ Filtration Rate Characterization

- Rosenmund Pocket Filtration


☐ Analytical Tools

- HPLC
- GC
- Karl Fischer
- Crystal 16




☐ Solid State Characterization Tools

- XRPD
- VTI (Vapor Sorption)
- Differential Scanning Calorimetry (DSC)
- Thermogravimetric Analysis (TGA)
- Raman Microscopy
- On-line Raman


In-situ microscopy (PVM)



Solubility and solute concentration (ATR-FTIR)






Chord length distribution (FBRM)




Polymorphic form (Raman)

- ❑ Preliminary assessment of a compound's propensity for growth
 - Quick regression of growth parameters gives preliminary read regarding whether crystallization kinetics are:
 - ❑ Fast posing potential challenges with need for rapid micro/meso-mixing
 - ❑ Slow creating potential time cycle challenges
- ❑ Rapid screening of batch crystallization conditions
 - Numerous alternative crystallization modes can be simulated in-silico, with only minimal lab-scale pre-investment required
- ❑ Optimization of target process parameters/recipes
 - Advanced process optimization algorithms enable more efficient process design and opportunities for increased process robustness
- ❑ Batch-to-Continuous Process Development
 - Lab-scale continuous crystallization capabilities are severely limited with existing lab infrastructure
 - Simulation offers an attractive tool for design optimization that is both time and resource efficient and can be based on existing batch process information



- Heavily seeded, isothermal crystallizations are executed to gain understanding of a compound's propensity for growth (vs. secondary nucleation for example)
- Comparison of growth rate to historical database enables rapid assessment of expected process cycle times and potential processing challenges (e.g. need for rapid mixing)



- Flow sheeting environment for building the crystallization
- Options for temperature, level and pressure controllers, as well as numerous sensors/ analyzers (ie. PSD)
- A convenient graphical user interface for parameter entry with numerous crystallization kernels built in

vessel (crystallizer_MSMPR)

Slurry level:

Temperature:

Kinetics information:

Growth and dissolution:

Specify


<input checked="" type="checkbox"/> Growth rate constant	<input checked="" type="radio"/> Uniform for entire array <input type="radio"/> Per element	m/s
	<input type="text" value="1.63985e-6"/>	
<input checked="" type="checkbox"/> Order with respect to supersaturation	<input checked="" type="radio"/> Uniform for entire array <input type="radio"/> Per element	
	<input type="text" value="1.17104"/>	
<input checked="" type="checkbox"/> Activation energy	<input checked="" type="radio"/> Uniform for entire array <input type="radio"/> Per element	J/mol
	<input type="text" value="16573"/>	
<input checked="" type="checkbox"/> Effective diffusivity correction factor	<input checked="" type="radio"/> Uniform for entire array <input type="radio"/> Per element	
	<input type="text" value="1"/>	

Secondary nucleation - attrition Growth and dissolution Agglomeration




Liquid phase reactions Initial conditions Initial PSD

Configuration Classification Primary nucleation Secondary nucleation - activated


OK Cancel Reset All Help



- ❑ Seeded isothermal crystallization experiments with a high seed loading in varying solvent mixtures
- ❑ Growth parameters estimated using gCrystal
 - gCrystal's graphical user interface and straightforward regression techniques simplifies growth parameter acquisition
- ❑ Understanding of growth kinetics helps guide additional crystallization development




$$G = k_1(S - 1)^{k_2}$$




0.4% seed @ 65C
Age at seed point
Staged cool down

Secondary Nucleation
then Growth



20% seed @ 70°C
Age at seed point
Staged cool down

Growth &
Agglomeration Only



5% seed @ 71°C
Staged cool down

Growth &
Agglomeration then
Secondary nucleation




0.4% seed @ 65°C

Age at seed point

Staged cool down

Secondary Nucleation
then Growth




20% seed @ 70°C

Age at seed point

Staged cool down




Growth &
Agglomeration Only



5% seed @ 71°C



Staged cool down

Growth &
Agglomeration then
Secondary nucleation





	D10	D50	D90
Model prediction - Growth & 2nd Nucleation	22.5	82.5	158.6
Experimental Data -	22.5	72.5	224.1

- Distribution broadening prevents accurate prediction of D90 results
- Broadening believed to be attributed to shape factors associated with laser diffraction techniques





	D10	D50	D90
Model prediction - Growth only	25.2	45.7	93.3
Model prediction - + Agglomeration	115.3	183.4	275
Experimental Data - No sonication	100.0	207.1	338.9

- SEM micrographs indicate minimal secondary nucleation
- Growth only model does not accurately predict the data
- Addition of agglomeration kernel improves PSD prediction significantly
- Degree of sonication during analysis impacts PSD significantly





	D10	D50	D90
Model prediction - Growth only	25.2	45.7	93.3
Model prediction - + Agglomeration	115.3	183.4	275
Experimental Data - No sonication	100.0	207.1	338.9

- SEM micrographs indicate minimal secondary nucleation
- Growth only model does not accurately predict the data
- Addition of agglomeration kernel improves PSD prediction significantly
- Degree of sonication during analysis impacts PSD significantly




	D10	D50	D90
Model prediction - Growth only	25.2	45.7	93.3
Model prediction - + Agglomeration	115.3	183.4	275
Experimental Data - No sonication	100.0	207.1	338.9

- SEM micrographs indicate minimal secondary nucleation
- Growth only model does not accurately predict the data
- Addition of agglomeration kernel improves PSD prediction significantly
- Degree of sonication during analysis impacts PSD significantly



	D10	D50	D90
Model prediction - Growth and 2nd Nucleation	32.3	58.1	104.8
Model prediction - + Agglomeration	64.9	129.2	206.5
Experimental Data -	29.2	118.0	264.3

- SEM micrographs indicate significant secondary nucleation, at the end of the crystallization
- Growth and 2nd nucleation only model does not accurately predict the data
- Model is not accurately capturing bimodal behavior



□ Agglomeration

- Agglomeration parameter regression was not functioning properly at the time of original investigation. Quantitative prediction of the agglomeration parameter should improve further agglomeration kernel performance.
- Impact of degree of sonication employed during analysis impairs agglomeration regression. Re-evaluation with a compound less prone to agglomeration is currently underway.

□ Shape effects in laser diffraction measurement

- Re-execute population balance model evaluation with more block like material
- Current investigations are planned towards repeating particle size measurement with alternative techniques less sensitive to particle shape (i.e. Coulter counter)
 - Techniques such as coulter counter are less frequently used within the industry and therefore do not likely serve to address the issue
 - Can data generated as output from gCrystal be treated so as to predict the broadening attributed to such shape effects? we shall see...

- Propensity for growth
 - Growth kinetics are rapidly and easily estimated using the gCrystal software environment and helpful in assessing a particular substance's crystallization tendencies.

- Particle size distribution prediction
 - Population balance models for prediction of particle size distribution can reasonably predict D50, provided the appropriate kernels are active and regressed, and alteration of the simulation parameters are rapid

 - Measurement error attributed to laser diffraction measurement techniques appears to be limiting the quantitative prediction from a particle size perspective, but qualitatively, model is well suited for determining mean size and active crystallization mechanisms.

