

The utility of the Artificial Stomach Duodenum system in *in vivo* oral absorption modeling

Kaoutar Abbou Oucherif

Advanced Modeling Forum, London , April 2016

Outline


- Study objectives
- Experimental setup
- gCOAS ASD in silico model
- Simulation results
- Oral absorption simulations
- Conclusions

Study Objectives

- Understand the extent of precipitation for different salts of a Lilly molecule in the GI tract
- Extract precipitation kinetic parameters from the ASD in vitro experiments to be used in the in vivo oral absorption model

Experimental setup

All experiments were done under dog fasted conditions



Nucleation Models

- **Classical nucleation:**

- $J_{prim} = \ln A_0 \left(\frac{-16\pi(\alpha\sigma)^3 v^2_0}{3k^3 T^3 \ln S^2} \right)$
- Parameters to be estimated: pre-exponential factor and the surface energy correction factor.

- **Power law kinetics:**

- $J_{prim} = \ln k_n \left(\frac{\Delta C}{\rho_c} \right)^n \exp\left(\frac{-E_{A,n}}{RT}\right)$
- Parameters to be estimated: nucleation coefficient, nucleation order, and the activation energy.

- **Custom kinetics:**

- E.g., extracting nucleation rates from probability distribution functions of induction time $P(t) = 1 - \exp(jV(t - t_g))$
- Input nucleation rate as a function of physiological parameters : $J_{prim} = \text{fcn}(\text{pH}, C_{bile})$

Growth model

- Surface integration limited growth:

$$\phi_{diss} = k_g (1 - S_{rel})^n$$


ϕ_{diss} : flux of an individual specie leaving the solid due to dissolution

k_g : surface integration rate constant

S_{rel} : relative supersaturation $\frac{C}{C^{sat}_{MW}}$

n : surface integration order

Dissolution



At the solid-liquid interface: $0 = \{\text{dissolution}\} + \{\text{reaction}\} + \{\text{diffusion}\}$

$$K_j = \prod_{i=1}^{NC} (C_i^S)^{v_{ij}}, \quad j = 1, \dots, NR$$

$$C_{i^*}^S = S_0$$

$$0 = R_{i^*} \delta_{i,i^*} + \sum_{j=1}^{NR} v_{ij} r_j - \frac{D_i}{h} (C_i^S - C_i^b), \quad i = 1, \dots, NC$$

gCOAS ASD in silico model

Global specifications

Please double click on the icon to enter:

1) Drug properties

Tablet

Please double click

1) Dosing schedule
2) Particle size dist

To change to a different
click on the icon and
MODEL"

global_specifications (Global_specifications)

Specification mode: Standard

Drug properties

Particle size grid setup

Specify

Solid type: Base

Dosed as: Freeform

☒ API name:

☒ Mass density: ☐ Uniform for entire array ☒ Per element

	Freeform	
Freeform	1.2	

g/cm3

☒ Molecular weight:

g/mol

☒ Stoichiometry: ☐ Uniform for entire array ☒ Per element

	Freeform				
API+	0				
API	1				
H2O	0				

-

☒ Intrinsic solubility: ☐ Uniform for entire array ☒ Per element

mg/ml


☒ pKa: ☐ Uniform for entire array ☒ Per element

	1st Dissociation	
1st Dissociation		


LOG10(#)

☒ Log P:

OK Cancel Reset all Help




Simulation setup




Experimental Results: concentration profiles

Gastric concentration profiles



Duodenal concentration profiles




Case studied: compound + citric acid


Lee Burns

Experimental Results: pH profiles

Stomach pH Profile




Duodenal pH Profile




Lee Burns

Experimental Results: pH profiles


Stomach pH Profile




Duodenal pH Profile



Estimating precipitation kinetics



Simulation results




Lack of fit test:


Weighted residual	χ^2 value (95%)
42.33	52.19

Model parameter	Final value	Initial guess	95% C.I	Standard deviation
Growth constant	0.55	0.36	1.42	0.70
Growth order	1.60	0.81	2.29	1.13
Nucleation coefficient	14.01	13.85	974.9	481.2
Nucleation order	2.61	2.77	399.3	197.1

Work flow of oral absorption modeling using gCOAS




A gPROMS-based Computational Oral Absorption Simulation Framework



K. Sugano, *Expert Opin. Drug Metab. Toxicol.* (2009) 5(3), pp. 259-293

Incorporating ASD results in gCOAS GI model



Final remarks

- The gCOAS ASD framework allows to extract quantitative precipitation kinetic parameters from ASD data
- The kinetic parameters for precipitation can be integrated in the gCOAS oral absorption model to better understand the impact of precipitation on bioavailability
- The ASD setup along with the gCOAS tools can further help investigate the impact of bile salts and pH on precipitation
- Custom precipitation kinetics can allow to account for changes in physiology (pH, bile concentration) on nucleation and growth